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The Selflab **EHR Ab** test is used for the detection of antibodies against ehrlichiosis in canine blood. An antibody (Ab) is a type of protein produced by the immune system in response to contact with an antigen.

**RECOMMENDATIONS AND PRECAUTIONS:**

1. Before using the product, read the instruction manual. Follow the instructions carefully. Otherwise, the results may be inaccurate.
2. Do not use the test kit beyond its expiration date.
3. Ensure all items in the test kit and specimens are at the same temperature before testing - this will take approximately 15-25 minutes.
4. The test should remain in the sealed pouch until use. Perform the test immediately after removing the test cassette from the pouch. Touch only the sides of the test cassette.
5. Do not use the test kit if the pouch is open or damaged.
6. Do not mix test kit components from different lot numbers.
7. This test kit is for external use only. Do not swallow any of the test kit components.
8. Avoid getting buffer solution in the eyes or on the skin. In case of accidental contact, rinse with plenty of water to avoid skin irritation. In case of eye contact, flush eyes with plenty of water for at least 15 minutes. If swallowed, drink plenty of water. Seek medical attention.
9. In the event of buffer solution, clean up spills thoroughly with a disinfectant.
10. Use only the buffer solution (from the tube provided in the kit) for testing. The use of other liquids (e.g., water) may yield inaccurate results.
11. Keep the test kit out of the reach of children.
12. The test kit is for single use only. Do not reuse any of the components of the test kit.
13. Dispose of the used test kit components in a sealed bag for mixed waste.
14. Wash your hands thoroughly before and after testing.

**Storage of Tests:**

Store the test kits in a temperature 2-30°C. Do not freeze, protect from direct sunlight.

**Sample Storage:**

The collected sample should be tested immediately.

**USE CAUTION:**

Samples are considered infectious material! Expiration date is indicated on the package. Tests are for animal use only.

**Performance:**


**Test specificity** is the ratio of the true negatives to the sum of the true negatives and false positives combined and was: for EHR 99.6% for this assay. It means that 99.6% of samples from healthy individuals with this test are marked as healthy / negative.

The **sensitivity** of the test is the ratio of the true positives to the sum of the true positives and the false negatives. The sensitivity: of EHR 99%, estimated in this test, means that such a percentage of sick individuals will be positively diagnosed.

**LIMITATIONS:**

1. The test result should not be the sole source of diagnosis, as with other diagnostic tests, the result should be interpreted by the veterinarian with the other clinical data.
2. False positive result may occur when:
  - there is a cross-reactivity, there is another infection in the organism,
  - the cured patient has persistent antibodies,
  - patient has an autoimmune disease.
3. False negative result may occur when:
  - the specimen has been sampled incorrectly or sample is too diluted,
  - antibodies levels in the sample are too low,
  - the sample is taken too early (serological window), the body has not yet had the time to produce antibodies.
4. Due to the possibility of false-positive or false-negative results, cassette test results may be confirmed by other methods.
5. Invalid test result may occur when:
  - the specimen was sampled incorrectly or sample is too diluted,
  - not enough drops of sample material have been placed in the test well "S".
6. If the symptoms persist despite the negative or invalid results, it is recommended to do the test again within 24-48h and contacting your veterinarian.
7. Positive test result may only indicate that the pet had contact with microorganisms and doesn't have to be equivalent of the disease development.
8. Negative result does not exclude a possibility of infection.
9. Positive result does not exclude a possibility of coexisting infection with different pathogens.
10. Each positive result should be consulted with a veterinarian.

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**EN Positive result:** there will be 2 lines, 1 in control field (C) and 1 in the test field (T). **Look closely! The presence of any line in test field (T), even very vague, means that the result is considered positive.**  
**Negative result:** only 1 clear line in control field (C), no line in test field (T).  
**Invalid result:** no color line in control field (C), the invalid result is not affected by the appearance of the test field (T) (presence or absence of a line in this field).

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**EN Test principle:**  
The test is performed by a lateral flow immunochromatographic method (to detect the presence or absence of the target Ag in the sample). The test cassette has a window (S) into which the sample is dropped and two fields: test field (T) and control field (C). When the sample is placed in the (S) window, the liquid will flow laterally into test field (T). If the antigen is present in the sample, a visible (T) line will appear.

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



### CONSTRUCTION OF ALL-IN-ONE CONTAINER

#### BEHÄLTER-BAU BUDOWA POJEMNIKA

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EN

Take all the items out of the box: gloves (put them on), lancet for puncture, disinfecting swab, razor, small adhesive patch, assay buffer container, all-in-one container (blue) - put them upright.

Choose the puncture spot (ear or paw). If necessary use the razor and shave the puncture spot.

Take the alcohol swab, take it out of the packaging and disinfect the puncture spot. Wait until the puncture spot will naturally dry.

Take the all-in-one container, twist the upper part anti-clockwise and take it off.

Take the container with the buffer, unscrew the cap and pour its contents into the all-in-one container.

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





or  
oder  
lub

EN

Twist off the lancet blue protective cap and pull it straight out.

Disinfect the insertion area. Press the lancet firmly into the insertion site and press the blue button until you hear a click.

Blood from the injection site should be drawn into into all-in-one container. If no blood appears, you can massage the area in the direction of the puncture.

After taking the sample, apply firm pressure to the puncture site. Apply a band-aid if needed.

Put the top of the container with the collected blood into the all-in-one container, close it and mix the blood well with the buffer.

Squeeze the all-in-one cap until you will hear a "click" sound to take it off the upper part.

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