SCHNELLDIAGNOSETESTS

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
- EN** **Positive result:** there will be 2 lines, 1 in control field (C) and 1 in the test field (T). **Look closely! The presence of any line in test field (T), even very vague, means that the result is considered positive.**
Negative result: only 1 clear line in control field (C), no line in test field (T).
Invalid result: no color line in control field (C), the invalid result is not affected by the appearance of the test field (T) (presence or absence of a line in this field).

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- EN** **Test principle:**
The test is performed by a lateral flow immunochromatographic method. The test cassette has a window (S) into which the sample is dropped and two fields: test field (T) and control field (C). When the sample is placed in the (S) window, the liquid will flow laterally into test field (T). If the antigen is present in the sample, a visible (T) line will appear.

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RAPID DIAGNOSTIC TESTS

The Selflab **CAV1-ICH Ag** test is used to detect adenovirus type 1 antigen in nasal secretions causing infectious canine hepatitis. Antigen (Ag) is a foreign substance that stimulates the body to produce antibodies.

RECOMMENDATIONS AND PRECAUTIONS:

- Before using the product, read the instruction manual. Follow the instructions carefully. Otherwise, the results may be inaccurate.
- Do not use the test kit beyond its expiration date.
- Ensure all items in the test kit and specimens are at the same temperature before testing - this will take approximately 15-25 minutes.
- The test should remain in the sealed pouch until use. Perform the test immediately after removing the test cassette from the pouch. Touch only the sides of the test cassette.
- Do not use the test kit if the pouch is open or damaged.
- Do not mix test kit components from different lot numbers.
- This test kit is for external use only. Do not swallow any of the test kit components.
- Avoid getting buffer solution in the eyes or on the skin. In case of accidental contact, rinse with plenty of water to avoid skin irritation. In case of eye contact, flush eyes with plenty of water for at least 15 minutes. If swallowed, drink plenty of water. Seek medical attention.
- In the event of buffer solution, clean up spills thoroughly with a disinfectant.
- Use only the buffer solution (from the tube provided in the kit) for testing. The use of other liquids (e.g., water) may yield inaccurate results.
- Keep the test kit out of the reach of children.
- The test kit is for single use only. Do not reuse any of the components of the test kit.
- Dispose of the used test kit components in a sealed bag for mixed waste.
- Wash your hands thoroughly before and after testing.

Storage of Tests:

Store the test kits in a temperature 2-30°C. Do not freeze, protect from direct sunlight.

Sample Storage:

The collected sample should be tested immediately.

USE CAUTION:

Samples are considered infectious material! Expiration date is indicated on the package. Tests are for animal use only.

Performance:

Test specificity is the ratio of the true negatives to the sum of the true negatives and false positives combined and was: for CAV1-ICH Ag 98.7% for this assay. It means that 98.7% of samples from healthy individuals with this test are marked as healthy / negative.

The sensitivity of the test is the ratio of the true positives to the sum of the true positives and the false negatives. The sensitivity: of CAV1-ICH Ag 93.6%, estimated in this test, means that such a percentage of sick individuals will be positively diagnosed.

LIMITATIONS:

- The test result should not be the sole source of diagnosis, as with other diagnostic tests, the result should be interpreted by the veterinarian with the other clinical data.
- False positive result may occur when:
 - there is a cross-reactivity, there is another infection in the organism.
- False negative result may occur when:
 - the specimen has been sampled incorrectly,
 - antigens levels in the sample are too low,
 - the sample has been taken at the wrong time - in the sampled material there's no antigen.
- Due to the possibility of false-positive or false-negative results, cassette test results may be confirmed by other methods.
- Invalid test result may occur when:
 - the specimen was sampled incorrectly,
 - not enough drops of sample material have been placed in the test well "S".
- If the symptoms persist despite the negative or invalid results, it is recommended to do the test again within 24-48h and contacting your veterinarian.
- A positive result may be only indicative of contact with the microorganism and doesn't have to be equivalent of the disease development.
- A negative result does not exclude a possibility of infection.
- Positive result does not exclude a possibility of coexisting infection with different pathogens.
- Each positive result should be consulted with a veterinarian.

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



EN

Take all the items out of the box: gloves (put them on), swab, assay buffer container (put them upright).

Place the supplied swab up to 1 cm into the dog's nose and collect a sample of secretions by rotating it several times along the mucous membrane inside the nasal cavity.

Twist off the assay tube cap.

Put the swab with the sample into the tube, next mix the sample with the buffer well.

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


EN

Take out the test cassette and place it on a flat and clean surface.

Using the pipette from the test pouch, take the sample.

Add 4 drops into the test's sample well (S).

Read the results after 5-10 minutes. **WARNING! Results read after 10 minutes may be invalid!**

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
EN

After use collect all the test components to an empty trash bag, close it precisely. Throw it out to mixed trash.

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selflab

Sigmed Sp. z o.o.
Mickiewicza 59,
47-253 Cisek
Poland

E-mail: selflab@selflab.info
Tel:
DE: + 49 301 20 74 561
PL/EN: +48 77 707 08 48
www.selflab.info